Oral-B®

Professional Care OxyJet® Center

OC 17 565 X

powered
by

BRAUN

Type 3719

Thank you for purchasing an Oral-B Professional Care OxyJet® Center.

If you have any questions, please call:

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
IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

Read all instructions before using.

DANGER

To reduce the risk of electrocution:

1. Do not place or store charger where it can fall or be pulled into a tub or sink.
2. Do not place charger in or drop into water or other liquid.
3. Do not reach for a charger that has fallen into water. Unplug immediately.
4. Do not use while bathing.

WARNING

To reduce the risk of burns, electrocution, fire, or injury to persons:

1. Close supervision is necessary when this product is used by, on, or near children.
2. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
3. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or if the charger has been dropped into water. Return the product to a service center for examination and repair.
4. Keep the cord away from heated surfaces.
5. Never drop or insert any object into any opening.

SAVE THESE INSTRUCTIONS

The Oral-B Professional Care OxyJet® Center has been developed with leading dental professionals to effectively remove plaque from all areas of your teeth and to attack plaque bacteria.

The Oral-B Professional Care toothbrush adds a gentle pulsating movement to its oscillation to create a deep-cleaning 3D brushing action. The toothbrush should be used twice a day for at least two minutes, the minimum brushing time recommended by dental professionals.

The Oral-B OxyJet® irrigator mixes air into water in the form of micro-bubbles. These micro-bubbles in the water stream attack plaque bacteria. At the same time, the irrigator removes remaining food particles from between your teeth and massages your gums. Always use the irrigator after brushing your teeth.

Please read the use instructions carefully before first use.

Important

The plug should always be accessible. Periodically check the cord for damage. A damaged or non-functioning appliance should no longer be used. Take it to a Braun Service Center.

Description

Toothbrush

- ① Charging base
- ② Charge indicator
- ③ Toothbrush handle (rechargeable)
- ④ On/off switch
- ⑤ Brushhead
- ⑥a Power Tip® brushhead
- ⑥ Brushhead compartment








OxyJet® irrigator

- ⑦ Jet compartment
- ⑧ Irrigator jet
- ⑨ Jet switch
- ⑩ Jet release button
- ⑪ On/off switch
- ⑫ Water pressure control dial
- ⑬ Handle switch (water stop and go)
- ⑭ Irrigator handle
- ⑮ Water container (20 fl.oz / 600 ml)

Specifications

Power supply: see bottom of housing

Power intake: 24 W

Using the toothbrush

Connecting and charging

The toothbrush is electrically safe and can be used in the bathroom without hesitation.

- Plug the appliance into an electrical outlet. Place the toothbrush handle ③ on the charging base ①.
- The charge indicator ② flashes while the handle is being charged.
- When the handle is fully charged, the charge indicator lights permanently. A full charge takes at least 16 hours and provides an operation time of approximately 45 minutes.
- For everyday use, the toothbrush handle can be stored on the charging base to maintain it at full power. Overcharging is impossible.

Battery maintenance

To maintain the maximum capacity of the rechargeable battery, the following procedure should be done every 6 months: Do not place the toothbrush handle back onto the charging base in order to discharge it by regular use. Once discharged, recharge the handle on the charging base.

Operating the on/off switch ④

The toothbrush features two speed options, the lower speed being ideal for sensitive areas in the mouth.

- Press the switch ④ to operate with normal speed.
- Press a second time to operate with reduced speed.
- Press a third time to switch the toothbrush off.


Removing plaque

- The toothbrush can be used with any leading brand of toothpaste.
- To avoid splashing, guide the brushhead to your teeth before switching on.
- Guide the brushhead slowly from tooth to tooth. Hold the brushhead in place for a few seconds before moving on the next tooth.
- Brush the gums as well as the teeth, first the outsides, then the insides and finally the chewing surfaces. Do not press too hard nor scrub, simply let the brush do all the work.

During the first days of using the toothbrush, your gums may bleed slightly. In general, bleeding should stop after a few days. Should it persist after two weeks, please consult your dentist or hygienist.

Professional timer

With a signal you can hear and feel, the timer invites you to brush quadrant by quadrant. After every 30 seconds, it signals by a very short stutter that it is time to move to the next quadrant. After two minutes - which is the recommended brushing time - the timer signals by a longer stutter.


Pressure sensor

For safety and optimum brushing results, the toothbrush has a pressure control system: If too much pressure is applied, the oscillating movement of the brushhead will continue, but the pulsation will stop. There is a noticeable change in feeling and sound when the pressure sensor is activated.

INDICATOR® bristles


The INDICATOR® bristles help you monitor and improve the brushing performance.

With proper brushing, twice a day for two minutes, the blue color on the INDICATOR® bristles will disappear halfway down within 3 months, signalling it is time to replace the brushhead.

Power Tip® brushhead

The Power Tip® brushhead is designed for cleaning between teeth as well as around bridges, crowns and implants. To use, position brushhead at 45° angle to gumline. Slowly guide the brushhead along the gumline, pausing between teeth for optimum results.


Replacement brushheads

Available at retail stores or Braun Service Centers.

After use

Rinse the brushhead thoroughly under running water for several seconds with the handle switched on.

Switch off the handle and remove the brushhead. Clean both parts separately under running water, then wipe them dry.


The environment





For environmental reasons, the toothbrush handle is provided with an «environmentally sound» nickel-hydrate rechargeable battery. The battery does not contain any harmful heavy metals. Nevertheless, in the interest of recycling, please do not dispose of the product in the household waste at the end of its useful life. Instead dispose of it at an authorized Braun Service Center or at an appropriate collection site.

However, if you prefer to dispose of the battery yourself, open the handle as shown, remove the battery and dispose of it according to local environmental regulations. Caution:

Opening the handle will definitely destroy it.

Note: You will not be able to place new batteries in the handle.

Using the OxyJet® irrigator

Before using for the first time, fill the water container and allow the water to run through completely by holding the irrigator handle over a washbasin.

1.
Attach one of the irrigator jets to the handle until it locks into place. Fill the water container with water and add mouthwash, if desired.

2.
Turn on the irrigator by pressing the on/off switch. (When doing so, the handle switch should be in the lower position = water stop).

3.
Select the water pressure with the pressure control dial (1 = gentle, 5 = strong). We recommend that you start with setting 1, especially if your gums tend to bleed or if children are using the irrigator.

4.
Bend over a washbasin and guide the jet to your teeth and gums. Slide the handle switch up to start the water stream. Keep your mouth slightly open to allow the water to run into the washbasin.




You can choose between two different jet functions by sliding the jet switch ⑨ down or up (before doing so, turn off the water stream by sliding the handle switch ⑬ down):

Rotating stream for general cleaning and gum massage (jet switch in lower position): The rotation of the water stream and its micro-bubbles is designed for easy penetration below the gumline.

Mono stream: straight water stream for site-specific cleaning (jet switch in upper position)

Max. operating time: 15 minutes
Cooling off time: 2 hours

To ensure the jet functions, a small sieve is incorporated in the lower end of the shaft. Should you notice a reduced water pressure, check the sieve for any deposits and remove them, e.g. by carefully tapping the shaft on a flat surface.

After use

Always empty the water container entirely to prevent bacteria from collecting in the standing water. If your dentist has recommended the use of a chlorhexidine solution with your irrigator, allow some water to run through the irrigator after use to prevent any deposit build-up.


Press the on/off switch to turn off the irrigator. It will automatically turn off after 10 minutes if you forget to turn it off manually. In this case, to turn it on again, you have to press the on/off switch twice.

To take off the irrigator jet, press the jet release button. The jets can be stored in the jet compartment. Wipe the handle dry and place it back onto its base.

Cleaning

Regularly clean the housing with a damp cloth. The water container and the jet compartment are removable and dishwasher-safe.

Subject to change without notice

FOR USA ONLY

Limited 2-Year Warranty

Braun warrants this appliance to be free of defects in material and workmanship for a period of two years from the date of original purchase at retail.

If the appliance exhibits such a defect, Braun will, at its option, repair or replace it without cost for parts and labor. The product must be carried in or shipped, prepaid and insured (recommended), to an authorized Braun Service Center. Proof of date of original purchase is required.

The cost to ship the product to an authorized Service Center is the responsibility of the consumer, the cost to return the product to the consumer is the responsibility of Braun.

This warranty does not cover finishes or normal wear, nor does it cover damage resulting from accident, misuse, dirt, water, tampering, unreasonable use, service performed or attempted by unauthorized service agencies, failure to provide reasonable and necessary maintenance, including recharging, or units that have been modified or used for commercial purposes.

ALL IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO TWO YEARS FROM DATE OF ORIGINAL PURCHASE. IN NO EVENT WILL BRAUN BE RESPONSIBLE FOR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE OF THIS PRODUCT.

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you.

This warranty gives you specific legal rights, and you may have other legal rights which vary from state to state.

How to obtain service in the USA

Should this product require service, please call 1-800-BRAUN 11 (272-8611) to be referred to the Authorized Braun Service Center closest to you.

To obtain service:

- A. Carry the product in to the authorized Braun Service Center of your choice or,
- B. Ship the product to the authorized Braun Service Center of your choice.
Pack the product well.
Ship the product prepaid and insured (recommended). Include a copy of your proof of purchase to verify warranty coverage, if applicable.

To obtain original Braun replacement parts and accessories:

Contact an authorized Braun Service Center. Please call 1-800-BRAUN 11 (272-8611) to be referred to the authorized Braun Service Center closest to you.

To order original Braun replacement parts and accessories (including replacement brushheads) by phone:

Please call 1-800-892-8534.

To speak to a Braun Consumer Service Representative:

Please call 1-800-984-9668.

For Canada only

Limited 2-Year Warranty

In the event this appliance fails to function within the specified warranty period because of defects in material or workmanship, and the consumer returns the unit to an authorized service centre, Braun Canada will, at its option either repair or replace the unit without additional charge to the consumer. This warranty does not cover any product which has been damaged by dropping, tampering, servicing performed or attempted by unauthorized service agencies, misuse or abuse. This warranty excludes units which have been modified or used for commercial purposes.

To the extent allowed by law, Braun Canada shall not be responsible for loss of the product, loss of time, inconvenience, commercial loss, special or consequential damages.

The provisions of this warranty are in addition to and not a modification of or subtraction from the statutory warranties and other rights and remedies contained in applicable provincial legislation.

How to obtain service in Canada

Should this product require service, please call 1-800-387-6657 to be referred to the authorized Braun Service Center closest to you.

To obtain service:

- A. Carry the product in to the authorized Braun Service Center of your choice or,
- B. Ship the product to the authorized Braun Service Center of your choice.
Pack the product well.
Ship the product prepaid and insured (recommended). Include a copy of your proof of purchase to verify warranty coverage, if applicable.

To obtain original Braun replacement parts & accessories:

Contact an authorized Braun Service Center. Please call 1-800-387-6657 to be referred to the Authorized Braun Service Center closest to you.

To speak to a Braun Service Representative:

Please call 1-800-387-6657.

24 Month Performance Guarantee

For USA only

Use the Oral-B Professional Care Center daily until you are fully familiar with this mechanical plaque removal instrument. If you are not convinced that the Oral-B Professional Care Center significantly reduces plaque, send it back for a full refund.

If not satisfied please send the Oral-B Professional Care Center within 24 months, post-paid and insured along with your name and address filled in below and the original sales slip indicating the purchase price and the date of purchase to:

24 Month Performance Guarantee

Braun Division,
The Gillette Company,
Sobin Park and A Street, MS 1F
South Boston, MA 02127

Sorry: No refunds will be made for any appliance damaged by accident, neglect or unreasonable use.

The Oral-B Professional Care Center must be shipped in its original packaging and be postmarked no later than 24 months after the date of purchase.

Please allow 6–8 weeks for handling.

Name of Regular Toothbrush

Your Name

Street

City

State

Zip Code

24 Month Performance Guarantee

For Canada only

Use the Oral-B Professional Care Center daily until you are fully familiar with this mechanical plaque removal instrument. If you are not convinced that the Oral-B Professional Care Center significantly reduces plaque, send it back for a full refund.

If not satisfied please send the Oral-B Professional Care Center within 24 months, post-paid and insured along with your name and address filled in below and the original sales slip indicating the purchase price and the date of purchase to:

24 Month Performance Guarantee

Braun Products
Gillette Canada Company
4 Robert Speck Parkway
Mississauga ON L4Z 4C5

Sorry: No refunds will be made for any appliance damaged by accident, neglect or unreasonable use.

The Oral-B Professional Care Center must be shipped in its original packaging and be postmarked no later than 24 months after the date of purchase.

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Gillette Manufactura, S.A. de C.V. /
Gillette Distribuidora, S.A. de C.V.